

Photosynthesis and Respiration

Plants make sugar, storing the energy of the sun into chemical energy, by the process of photosynthesis. When they require energy, they can tap the stored energy in sugar by a process called cellular respiration.

The process of photosynthesis involves the use of light energy to convert carbon dioxide and water into sugar, oxygen, and other organic compounds. This process is often summarized by the following reaction:



Cellular respiration refers to the process of converting the chemical energy of organic molecules into a form immediately usable by organisms. Glucose may be oxidized completely if sufficient oxygen is available by the following equation:



All organisms, including plants and animals, oxidize glucose for energy. Often, this energy is used to convert ADP and phosphate into ATP.

OBJECTIVES

In this experiment, you will

- use an O₂ Gas Sensor to measure the amount of oxygen gas consumed or produced by a plant during respiration and photosynthesis.
- use a CO₂ Gas Sensor to measure the amount of carbon dioxide consumed or produced by a plant during respiration and photosynthesis.
- determine the rate of respiration and photosynthesis of a plant.




Figure 1

MATERIALS

| | |
|-------------------------------------|-----------------------------|
| LabPro or CBL 2 interface | 250-mL respiration chamber |
| TI Graphing Calculator | plant leaves |
| DataMate program | 500-mL tissue culture flask |
| Vernier O ₂ Gas Sensor | lamp |
| Vernier CO ₂ Gas Sensor | aluminum foil |
| CO ₂ -O ₂ Tee | forceps |

PROCEDURE

1. Plug the O₂ Gas Sensor into Channel 1 and the CO₂ Gas Sensor into Channel 2 of the LabPro or CBL 2 interface. Use the link cable to connect the TI Graphing Calculator to the interface. Firmly press in the cable ends.
2. Turn on the calculator and start the DATAMATE program. Press to reset the program.
3. Set up the calculator and interface for an O₂ Gas Sensor and CO₂ Gas Sensor.
 - a. Select SETUP from the main screen.
 - b. If the calculator displays an O₂ Gas Sensor in CH 1 and a CO₂ Gas Sensor in CH2, proceed directly to Step 4. If it does not, continue with this step to set up your sensors manually.
 - c. Press to select CH 1.
 - d. Select OXYGEN GAS from the SELECT SENSOR menu.
 - e. Select parts per thousand (PPT) as the unit.
 - f. Press once, then press to select CH2.
 - g. Select CO2 GAS from the SELECT SENSOR menu.
 - h. Select parts per thousand (PPT) as the unit.
4. Set up the data-collection mode.
 - a. To select MODE, press (the up arrow key) twice and press .
 - b. Select TIME GRAPH from the SELECT MODE menu.
 - c. Select CHANGE TIME SETTINGS from the TIME GRAPH SETTINGS menu.
 - d. Enter "15" as the time between samples in seconds.
 - e. Enter "40" as the number of samples (data will be collected for 10 minutes).
 - f. Select OK twice to return to the main screen.
5. Obtain several leaves from the resource table and blot them dry, if damp, between two pieces of paper towel.
6. Place the leaves into the respiration chamber, using forceps if necessary. Wrap the respiration chamber in aluminum foil so that no light reaches the leaves.
7. Insert the CO₂-O₂ Tee into the neck of the respiration chamber. Place the O₂ Gas Sensor into the CO₂-O₂ Tee as shown in Figure 1. Insert the sensor snugly into the Tee. The O₂ Gas Sensor should remain vertical throughout the experiment. Place the CO₂ Gas Sensor into the Tee directly across from the respiration chamber as shown in Figure 1. Gently twist the stopper on the shaft of the CO₂ Gas Sensor into the chamber opening. Do not twist the shaft of the CO₂ Gas Sensor or you may damage it.

8. Wait two minutes, then select START to begin data collection. Data will be collected for 10 minutes.
9. When data collection has finished, remove the aluminum foil from around the respiration chamber.
10. Fill the tissue culture flask with water and place it between the lamp and the respiration chamber. The flask will act as a heat shield to protect the plant leaves.
11. Turn the lamp on. Place the lamp as close to the leaves as reasonable. Do not let the lamp touch the tissue culture flask.
12. Press to view the graph of O₂ GAS VS. TIME. Sketch a copy of your graph in the Graph section below. When finished, press to return to the graph menu.
Press once, then press to view the graph of CO₂ GAS VS. TIME. Sketch a copy of your graph in the Graph section below. When finished, press to return to the graph menu. Select MAIN SCREEN from the graph menu.
13. Perform a linear regression to calculate the rate of respiration/photosynthesis.
 - a. Select ANALYZE from the main screen.
 - b. Select CURVE FIT from the ANALYZE OPTIONS menu.
 - c. Select LINEAR (CH 1 VS TIME) from the CURVE FIT menu.
 - d. The linear-regression statistics for these two lists are displayed for the equation in the form:
$$Y=A*X+B$$
 - e. Enter the value of the slope, A, as the rate of respiration/photosynthesis in Table 1.
 - f. Press to view a graph of the data and the regression line.
 - g. Press to return to the ANALYZE menu.
 - h. Repeat Steps 13b – 13g to calculate the respiration/photosynthesis rate using the data from the CO₂ Gas Sensor (CH 2 VS TIME).
 - i. Select RETURN TO MAIN SCREEN from the ANALYZE menu.
14. Repeat Steps 8 – 13 to collect data with the plant exposed to light.
15. Remove the plant leaves from the respiration chamber, using forceps if necessary. Clean and dry the respiration chamber.

DATA


| Table 1 | | |
|--------------|---|--|
| Leaves | O ₂ rate of production/consumption (ppt/s) | CO ₂ rate of production/consumption (ppt/s) |
| In the dark | | |
| In the light | | |

GRAPHS

Darkness



O₂ Gas vs. Time



CO₂ Gas vs. Time

Light



O₂ Gas vs. Time



CO₂ Gas vs. Time

QUESTIONS


1. Were either of the rate values for CO₂ a positive number? If so, what is the biological significance of this?
2. Were either of the rate values for O₂ a negative number? If so, what is the biological significance of this?
3. Do you have evidence that cellular respiration occurred in leaves? Explain.
4. Do you have evidence that photosynthesis occurred in leaves? Explain.
5. List five factors that might influence the rate of oxygen production or consumption in leaves. Explain how you think each will affect the rate?

EXTENSIONS

1. Design and perform an experiment to test one of the factors that might influence the rate of oxygen production or consumption in Question 5.
2. Compare the rates of photosynthesis and respiration among various types of plants.

TEACHER INFORMATION**Photosynthesis and Respiration**

1. Spinach leaves purchased from a grocery store work very well and are readily available any time of the year. For best results, keep the leaves cool until they are to be used. Just before use, expose the leaves to bright light for 5 minutes.
2. A fluorescent ring lamp works very well since it bathes the plant in light from all sides and it gives off very little heat. When using a ring lamp as shown below, it is not necessary to use a heat shield.




3. If tissue culture flasks are not available, a beaker or flask of water will also work. The tissue culture flask is very thin, however, and will allow leaves to receive much more light from the same lamp.
4. To extend the life of the O₂ Gas Sensor, always store the sensor upright in the box in which it was shipped.
5. The waiting time before taking data may need to be adjusted depending on the rate of diffusion of the oxygen gas and the carbon dioxide gas. Monitor the gas concentrations and start collecting data when the levels of gas begin to move in the correct direction.
6. The stopper included with the CO₂ Gas Sensor is slit to allow easy application and removal from the probe. When students are placing the probe in the CO₂-O₂ Tee, they should gently twist the stopper into the adapter opening. Warn the students not to twist the probe shaft or they may damage the sensing unit.
7. To conserve battery power, we suggest that AC Adapters be used to power the interfaces rather than batteries when working with the CO₂ Gas Sensor. An AC Adapter is shipped with each LabPro interface at the time of purchase. If you are using the CBL 2, you can purchase a Vernier AC Adapter for \$10 (order code-IPS).

SAMPLE RESULTS


| Table 1 | | |
|--------------|---|--|
| Leaves | O ₂ rate of production/consumption (ppt/s) | CO ₂ rate of production/consumption (ppt/s) |
| In the dark | -0.0023 | 0.00065 |
| In the light | 0.0045 | -0.00126 |

GRAPHS

Darkness




O₂ Gas vs. Time




CO₂ Gas vs. Time

Light



O₂ Gas vs. Time



CO₂ Gas vs. Time

ANSWERS TO QUESTIONS

1. The CO₂ rate value for leaves in the dark was a positive number. The biological significance of this is that CO₂ is produced during respiration. This causes the concentration of CO₂ to increase, as sugar is oxidized and broken into CO₂, water, and energy.
2. The O₂ rate value for leaves in the dark was a negative number. The biological significance of this is that O₂ is consumed during cellular respiration. This causes the concentration of O₂ to decrease as glucose is oxidized for energy.
3. Yes, cellular respiration occurred in leaves, since O₂ decreased when leaves were in the dark and photosynthesis was not possible.

4. Yes, photosynthesis occurred in leaves, since O₂ increased when leaves were exposed to light.
5. Answers may vary. They might include:
 - a. A greater number of leaves should increase the rate, since there are more chloroplasts to undergo photosynthesis and more cells to require energy through cellular respiration.
 - b. A greater light intensity will increase the rate of photosynthesis. It may not affect the rate of cellular respiration, however.
 - c. A cooler room may decrease both rates, as cellular metabolism decreases in cooler weather.
 - d. Facing the top of the leaves toward the light should increase the rate of photosynthesis, since the chloroplasts are closer to the light source.
 - e. If the plants overheat due to the heat from the lamp, they may wilt and stop functioning. This will decrease all rates.
 - f. If there are too many leaves, diffusion may be restricted and prevent accurate readings. This may apparently decrease both rates.