


Acid-Base Titration

A titration is a process used to determine the volume of a solution needed to react with a given amount of another substance. In this experiment, you will titrate hydrochloric acid solution, HCl, with a basic sodium hydroxide solution, NaOH. The concentration of the NaOH solution is given and you will determine the unknown concentration of the HCl. Hydrogen ions from the HCl react with hydroxide ions from the NaOH in a one-to-one ratio to produce water in the overall reaction:



When an HCl solution is titrated with an NaOH solution, the pH of the acidic solution is initially low. As base is added, the change in pH is quite gradual until close to the equivalence point, when equimolar amounts of acid and base have been mixed. Near the equivalence point, the pH increases very rapidly, as shown in Figure 1. The change in pH then becomes more gradual again, before leveling off with the addition of excess base.

In this experiment, you will use a pH Sensor to monitor pH as you titrate. The region of most rapid pH change will then be used to determine the equivalence point. The volume of NaOH titrant used at the equivalence point will be used to determine the molarity of the HCl.




Figure 1

OBJECTIVES

In this experiment, you will

- Use a pH Sensor to monitor changes in pH as sodium hydroxide solution is added to a hydrochloric acid solution.
- Plot a graph of pH vs. volume of sodium hydroxide solution added.
- Use the graph to determine the equivalence point of the titration.
- Use the results to calculate the concentration of the hydrochloric acid solution.

MATERIALS

Materials for *both* Method 1 (buret) and Method 2 (Drop Counter)

TI-Nspire handheld or computer and TI-Nspire software	250 mL beaker
data-collection interface	magnetic stirrer (if available)
Vernier pH Sensor	stirring bar or Microstirrer (if available)
HCl solution, unknown concentration	wash bottle
~0.1 M NaOH solution	distilled water
pipet bulb or pump	ring stand
	utility clamp

Materials required *only* for Method 1 (buret)

50 mL buret	2nd utility clamp
10 mL pipet	2nd 250 mL beaker

Materials required *only* for Method 2 (Drop Counter)

Vernier Drop Counter	100 mL beaker
60 mL reagent reservoir	10 mL graduated cylinder
5 mL pipet or graduated 10 mL pipet	

CHOOSING A METHOD

Method 1 has the student deliver volumes of NaOH titrant from a buret. After titrant is added, and pH values have stabilized, the student is prompted to enter the buret reading manually and a pH-volume data pair is stored

Method 2 uses a Vernier Drop Counter to take volume readings. NaOH titrant is delivered drop by drop from the reagent reservoir through the Drop Counter slot. After the drop reacts with the reagent in the beaker, the volume of the drop is calculated, and a pH-volume data pair is stored.

METHOD 1: MEASURING VOLUME USING A BURET

1. Obtain and wear goggles.
2. Use a pipet bulb (or pipet pump) to pipet 10 mL of the HCl solution into a 250 mL beaker. Add 50 mL of distilled water. **CAUTION:** *Handle the hydrochloric acid with care. It can cause painful burns if it comes in contact with the skin.*
3. Place the beaker on a magnetic stirrer and add a stirring bar. If no magnetic stirrer is available, you need to stir with a stirring rod during the titration.
4. Connect the pH Sensor to the data-collection interface. Connect the interface to the TI-Nspire handheld or computer.
5. Set up the data-collection mode.
 - a. Choose New Experiment from the  Experiment menu.
 - a. Choose Collection Mode ► Events with Entry from the  Experiment menu.
 - b. Enter **Volume** as the Name and **mL** as the Units. Select OK




Figure 2

6. Use a utility clamp to suspend a pH Sensor on a ring stand as shown in Figure 2. Position the pH Sensor in the HCl solution and adjust its position so that it is not struck by the stirring bar.
7. Obtain a 50 mL buret and rinse the buret with a few mL of the ~ 0.1 M NaOH solution. Dispose of the rinse solution as directed by your teacher. Use a utility clamp to attach the buret to the ring stand as shown in Figure 2. Fill the buret a little above the 0.00 mL level of the buret with ~ 0.1 M NaOH solution. Drain a small amount of NaOH solution so it fills the buret tip *and* leaves the NaOH at the 0.00 mL level of the buret. Record the precise concentration of the NaOH solution in your data table. **CAUTION:** *Sodium hydroxide solution is caustic. Avoid spilling it on your skin or clothing.*
8. You are now ready to perform the titration. This process is faster if one person manipulates and reads the buret while another person operates and enters volumes.
 - a. Click the Start button (▶) to prepare to collect data.
 - b. Before you have added any drops of NaOH solution, click the Keep button (📷) and enter 0 as the buret volume in mL. Select OK to store the first data pair for this experiment.
 - c. Add the next increment of NaOH titrant (enough to raise the pH about 0.15 units). When the pH stabilizes, click the Keep button (📷), enter the current buret reading (to the nearest 0.01 mL), and then select OK. You have now saved the second data pair for the experiment.
 - d. Continue adding NaOH solution in increments that raise the pH by about 0.15 units and enter the buret reading after each increment. When a pH value of approximately 3.5 is reached, change to a one-drop increment. Enter a new buret reading after each increment. **Note:** It is important that all increment volumes in this part of the titration be equal; that is, one-drop increments.

- e. After a pH value of approximately 10 is reached, again add larger increments that raise the pH by about 0.15 pH units, and enter the buret level after each increment.
 - f. Continue adding NaOH solution until the pH value remains constant.
9. Stop data collection () when you have finished collecting data.
 10. Examine the data on the displayed graph of pH vs. volume to find the *equivalence point*—that is the largest increase in pH upon the addition of 1 drop of NaOH solution. To examine the data pairs on the displayed graph, click any point on the graph. Use  and  to trace the data. Move to the region of the graph with the largest increase in pH. Find the NaOH volume just *before* this jump. Record this value in the data table. Then record the NaOH volume *after* the drop producing the largest pH increase was added. **Note:** Another method for determining the equivalence-point volume is described in the Alternate Equivalence Point Method of this experiment.
 11. (optional) Print a copy of the graph of pH vs. volume.
 12. Dispose of the beaker contents as directed by your teacher. Rinse the pH Sensor and return it to the pH storage solution.

METHOD 2: MEASURING VOLUME WITH A DROP COUNTER

1. Obtain and wear goggles.
2. Add 40 mL of distilled water to a 100 mL beaker. Use a pipet bulb (or pipet pump) to pipet 5.00 mL of the HCl solution into the 100 mL beaker with distilled water. **CAUTION:** *Handle the hydrochloric acid with care. It can cause painful burns if it comes in contact with the skin.*
3. Obtain approximately 40 mL of ~0.1 M NaOH solution in a 250 mL beaker. Record the precise NaOH concentration in your data table. **CAUTION:** *Sodium hydroxide solution is caustic. Avoid spilling it on your skin or clothing.*
4. Obtain the plastic 60 mL reagent reservoir. **Note:** The bottom valve will be used to open or close the reservoir, while the top valve will be used to finely adjust the flow rate. For now, close both valves by turning the handles to a horizontal position.

Rinse it with a few mL of the ~0.1 M NaOH solution. Use a utility clamp to attach the reagent reservoir to the ring stand. Add the remainder of the NaOH solution to the reagent reservoir.

Drain a small amount of NaOH solution into the 250 mL beaker so it fills the reservoir's tip. To do this, turn both valve handles to the vertical position for a moment, then turn them both back to horizontal.
5. Connect the pH Sensor to the data-collection interface. Lower the Drop Counter onto a ring stand and connect it to the data-collection interface. Connect the interface to the TI-Nspire handheld or computer.
6. Calibrate the Drop Counter so that a precise volume of titrant is recorded in units of milliliters.
 - a. Choose Set Up Sensors ► Calibrate ► Drop Counter from the  Experiment menu:

- If you have previously calibrated the drop size of your reagent reservoir and want to continue with the same drop size, select Manual Entry. Enter the value for Drops/mL, then select OK. Proceed directly to Step 7.
 - If you want to perform a new calibration, choose Single Point.
 - b. Place a 10 mL graduated cylinder directly below the slot on the Drop Counter, lining it up with the tip of the reagent reservoir.
 - c. Open the bottom valve on the reagent reservoir (vertical). Keep the top valve closed (horizontal).
 - d. Slowly open the top valve of the reagent reservoir so that drops are released at a slow rate (~1 drop every two seconds). You should see the drops being counted on the screen.
 - e. When the volume of NaOH solution in the graduated cylinder is between 9 and 10 mL, close the bottom valve of the reagent reservoir.
 - f. Enter the precise Volume of NaOH. Select Keep. Record the number of drops/mL displayed on the screen for possible future use. Select OK.
7. Discard the NaOH solution in the graduated cylinder as indicated by your instructor and set the graduated cylinder aside.
8. Assemble the apparatus.
 - a. Place the magnetic stirrer on the base of the ring stand.
 - b. Insert the pH Sensor through the large hole in the Drop Counter.
 - c. Attach the Microstirrer to the bottom of the pH Sensor, as shown in the small picture. Rotate the paddle wheel of the Microstirrer and make sure that it does not touch the bulb of the pH Sensor.
 - d. Adjust the positions of the Drop Counter and reagent reservoir so they are both lined up with the center of the magnetic stirrer.
 - e. Lift up the pH Sensor, and slide the beaker containing the HCl solution onto the magnetic stirrer. Lower the pH Sensor into the beaker.
 - f. Adjust the position of the Drop Counter so that the Microstirrer on the pH Sensor is just touching the bottom of the beaker.
 - g. Adjust the reagent reservoir so its tip is just above the Drop Counter slot.
9. Turn on the magnetic stirrer so that the Microstirrer is stirring at a fast rate.
10. You are now ready to begin collecting data. Check to see that the pH value is between 1.5 and 2.5. Start data collection (▶). No data will be collected until the first drop goes through the Drop Counter slot. Fully open the bottom valve—the top valve should still be adjusted so drops are released at a rate of about 1 drop every 2 seconds. When the first drop passes through the Drop Counter slot, check the graph to see that the first data pair was recorded.
11. Continue watching your graph to see when a large increase in pH takes place—this will be the equivalence point of the reaction. When this jump in pH occurs, let the titration proceed for several more milliliters of titrant, then stop data collection (◻) to view a graph of pH vs. volume. Turn the bottom valve of the reagent reservoir to a closed (horizontal) position.
12. Dispose of the beaker contents as directed by your teacher.



13. Examine the data on the displayed graph of pH vs. volume to find the *equivalence point*. To examine the data pairs on the displayed graph, click any point on the graph. Use ► and ◀ to trace the data. Move to the region of the graph with the largest increase in pH. Find the NaOH volume just *before* this jump. Record this value in the data table. Then record the NaOH volume *after* the drop producing the largest pH increase was added. **Note:** Another method for determining the equivalence-point volume is described in the Alternate Equivalence Point Method of this experiment.
14. (optional) Print copies of the graph.
15. If time permits, repeat the procedure. Be sure to Store the Latest Run (☒) before starting a new collection.

ALTERNATE EQUIVALENCE POINT METHOD

An alternate way of determining the precise equivalence point of the titration is to take the first and second derivatives of the pH-volume data.

1. Determine the peak value on the first derivative vs. volume plot.
 - a. Choose New Calculated Column from the ☒ Data menu.
 - b. Enter **1st Deriv** as the Name, **D1** as the Short Name, and leave the Units field blank.
 - c. Enter **derivative(pH, Volume, 1, 1)** as the Expression. Select OK.
 - d. Choose Select Y-axis Columns ► 1st Deriv from the ↙ Graph menu.
 - e. Examine the graph to determine the volume at the peak value of the first derivative.
2. Determine the zero value on the second derivative vs. volume plot.
 - a. Choose New Calculated Column from the ☒ Data menu.
 - b. Enter **2nd Deriv** as the Name, **D2** as the Short Name, and leave the Units field blank.
 - c. Enter **derivative(pH, Volume, 2, 1)** as the Expression. Select OK.
 - d. Choose Select Y-axis Columns ► 2nd Deriv from the ↙ Graph menu.
 - e. Examine the graph to determine the volume when the second derivative equals approximately zero. Note: You may have to estimate this value as it typically does not occur on an exact data point.

PROCESSING THE DATA

1. Use your printed graph and data table to confirm the volume of NaOH titrant you recorded *before* and *after* the largest increase in pH values upon the addition of 1 drop of NaOH solution.
2. Determine the volume of NaOH added at the equivalence point. To do this, add the two NaOH values determined above and divide by two.
3. Calculate the number of moles of NaOH used.
4. See the equation for the neutralization reaction given in the introduction. Determine the number of moles of HCl used.
5. Recall that you pipeted out 10.0 mL of the unknown HCl solution for the titration. Calculate the HCl concentration.

DATA AND CALCULATIONS

	Trial 1	Trial 2
Concentration of NaOH	M	M
NaOH volume added before largest pH increase	mL	mL
NaOH volume added after largest pH increase	mL	mL
Volume of NaOH added at equivalence point	mL	mL
Moles NaOH	mol	mol
Moles HCl	mol	mol
Concentration of HCl	mol/L	mol/L
Average [HCl]		M