



## Datalogging with Handheld Technology.

**KS3 Science**  
**Unit 7C**  
**Environment and feeding relationships**  
**Measuring physical environmental factors**

This year 7 unit provides many opportunities for datalogging with a variety of sensors, as can be seen from the following extract in the QCA specification:

How do environments vary?			
<i>Learning Objectives</i>	<i>Possible Teaching Activities</i>	<i>Learning Outcomes Pupils</i>	<i>Points To Note</i>
<ul style="list-style-type: none"> <li>♦ that some animals are adapted to daily changes in their habitat</li> <li>♦ how to measure and record changes in environmental factors</li> <li>♦ how to interpret patterns in data</li> </ul>	<ul style="list-style-type: none"> <li>♦ Ask pupils to predict how physical environmental factors around the school, <i>eg light intensity, temperature, humidity, noise levels</i>, would change over a 24-hour period and how they could measure the changes. With pupils, set up instruments, <i>eg datalogging equipment with a light probe, automatic weather station, temperature and sound sensors</i>, to monitor changes. Provide pupils with data about environmental changes around the school over a 24-hour period and help them to describe what these show and to identify links between the different changes.</li> </ul>	<ul style="list-style-type: none"> <li>♦ describe changes in physical environmental factors, <i>eg temperature, light intensity</i>, over a 24-hour period</li> <li>♦ interpret data about daily changes and explain in simple terms, <i>eg beginning to get dark</i></li> <li>♦ relate changes in variables, <i>eg light and temperature</i>, to each other</li> <li>♦ relate plant and animal activity to environmental changes</li> </ul>	<ul style="list-style-type: none"> <li>♦ There is an opportunity to use datalogging equipment. Secondary sources, <i>eg newspapers, Ceefax and geography weather stations</i>, could be used to confirm readings.</li> <li>♦ Schools near the coast and by tidal rivers have an opportunity to focus on adaptations to changes in these habitats.</li> <li>♦ Schools with CCTV may be able to monitor animal activity at night.</li> </ul>

### **Safety**

- ♦ teachers will need to check pupils' plans for health and safety before practical work begins.
- ♦ all off-site visits must be carried out in accordance with school/LEA guidelines.

### **Datalogging kit**

- ♦ TI-73, TI-82, TI-83 or TI-83plus graphing calculator.
- ♦ CBL or CBL2, with **fresh alkaline cells**.
- ♦ temperature, light, relative humidity sensors: microphone, (pH sensor, dissolved oxygen sensor).
- ♦ TI-GRAPH LINK™ software and cable to link to computer.
- ♦ chembio application (TI-83 plus) or program group (TI-73, TI-82 and TI-83).
- ♦ TI InterActive!™ or spreadsheet application.

### **Apparatus required**


- ♦ plastic bag/waterproof housing for CBL/CBL2.
- ♦ a variety of stands, clamps, bluetac etc. to keep sensors in position.
- ♦ flags/markers for seashore and riverside work.

### **Useful web sites**

[www.ti.com/calc/docs/graph.htm](http://www.ti.com/calc/docs/graph.htm)  
[www.vernier.com](http://www.vernier.com)  
[www.oxford-educational.co.uk](http://www.oxford-educational.co.uk)  
[www.qca.org.uk](http://www.qca.org.uk)

**Activity:****Measuring physical environmental factors.**

Because of the great variety of possible situations and habitats, this section will address the use of the CBL/CBL2 with multiple sensors over an extended period of time and will not concentrate on any specific activity. A combination of temperature, light and relative humidity sensors is appropriate for a school based or other terrestrial habitat. For small bodies of freshwater or rock pools, a pH sensor can be used in combination with temperature and light sensors. It is probably not appropriate to include a microphone in a multiple sensor set up because of the relatively long period of time between samples. Noise/sound levels can be investigated separately using a regime with more frequent sampling points.



In this example, lists containing captured data (time, **temperature**, light and **relative humidity**) have been transferred from the calculator to TI InterActive!™ and used to generate a graph. Lists can also be converted into .txt files using TI-GRAPH LINK™ and imported into spreadsheet applications.

**Technique**

1. Attach the calculator to the CBL/CBL2.
2. Open the chembio application/program and press enter to get to the following screen:-

```

***MAIN MENU***
1:SET UP PROBES
2:COLLECT DATA
3:VIEW GRAPH
4:VIEW DATA
5:FIT CURVE
6:RETRIEVE DATA
7:QUIT

```

3. Select SET UP PROBES and follow the on screen prompts to place temperature sensor in channel 1, light sensor in channel 2 and relative humidity sensor in channel 3.
4. On the CALIBRATION screen, select 1: USE STORED.
5. Select 2: COLLECT DATA from the MAIN MENU.
6. Select 2: TIME GRAPH.

7. For a 24 hour sampling period, sampling every 15 minutes or 30 minutes works well. For 15 minutes enter 900 or  $(15 \times 60)$  seconds. To record a data point at hour 24 (i.e. 86400 seconds), enter 97 or  $(24 \times 4) + 1$  samples; this will actually show an EXPERIMENT LENGTH of 87300.0S on the screen.
8. Select 1: USE TIME SETUP.
9. PRESS [ENTER] TO BEGIN COLLECTING DATA.

The calculator will now display the following screen:-


```

PERFORMING
EXPERIMENT...
WHEN CBL SHOWS
DONE, CHOOSE
RETRIEVE DATA
ON NEXT SCREEN.
[ENTER]

```

*The CBL screen can be used to trace the progress of the experiment. With the CBL2, however, it is necessary to keep a separate record of time as there is no screen.*

10. Detach the calculator, leaving only the CBL/CBL2 and sensors on site.
11. Press ENTER and then 7: QUIT. At this stage the TI-83 plus will show an ERR: INVALID DIM message: however, just press enter again to quit.
12. After sampling has finished ( **and it is important for the CBL/CBL2 to have finished collecting data**), re-attach the calculator, enter the chembio application/program and select 6: RETRIEVE DATA from the MAIN MENU. After the MAKE SURE CBL SHOWS DONE message, the screen will show TIME IN L1, DATA IN L2, L3, L4.



Pressing enter repeatedly will show time v. temperature, time v. light intensity and time v. relative humidity graphs in turn. Using TI-GRAPH LINK™, the screens can be 'grabbed' and either printed or copied to the clipboard. The lists of data can also be transferred to TI InterActive!™ or saved as .txt files and then opened in spreadsheet applications.

### Possible Extension Work.

- ◆ Use of other sensors (pH and dissolved oxygen) in freshwater studies.
- ◆ Investigation of noise levels with frequent sampling points.

### Links

- ◆ There are further opportunities to measure physical environmental factors in *Unit 8D Ecological relationships*.
- ◆ An investigation into noise levels forms an important part of *ICT Unit 7 Measuring physical data*.