

Energy Content of Foods

You are a lab technician working for NASA. Recently you were given the job of deciding what type of foods should be included in the next space mission. Four food types have been selected as possible snacks for the astronauts. You need to determine which of these four food choices has the highest energy content while adding the least amount of mass to the mission.

Your team will test two of the food types using a method known as calorimetry. During this process, you will burn a food sample positioned below a can containing a given amount of cold water. The water temperature will be monitored during the experiment using a temperature probe. By calculating the temperature change of the water, you will determine how much energy was released when the food sample burned.




FIGURE 1

OBJECTIVES

- **Measure** temperature changes.
- **Calculate** energy changes using specific heat.
- **Infer** the energy content of food.
- **Relate** energy content to types of food.
- **Evaluate** whether the nutrition labels are accurate.

MATERIALS

- can, small
- food samples (2)
- matches
- water, cold
- wooden splint

Energy Content of Foods *continued*

EQUIPMENT

- food holder (see **Figure 1**)
- graduated cylinder, 100 mL
- LabPro or CBL2 interface
- stirring rods (2)
- ring stand and 4-inch ring
- TI graphing calculator
- utility clamp and slit stopper
- Vernier temperature probe

SAFETY

- Wear safety goggles when working around chemicals, acids, bases, flames, or heating devices. Contents under pressure may become projectiles and cause serious injury.
- Secure loose clothing, and remove dangling jewelry. Do not wear open-toed shoes or sandals in the lab.
- Wear an apron or lab coat to protect your clothing when working with chemicals.
- In order to avoid burns, wear heat-resistant gloves whenever instructed to do so.
- If you are unsure of whether an object is hot, do not touch it.
- Avoid wearing hair spray or hair gel on lab days.
- Whenever possible, use an electric hot plate as a heat source instead of an open flame.
- Never return unused chemicals to the original container; follow instructions for proper disposal.

Procedure**EQUIPMENT PREPARATION**

1. Obtain and wear goggles.
2. Plug the temperature probe into Channel 1 of the LabPro or CBL 2 interface. Use the link cable to connect the TI graphing calculator to the interface. Firmly press in the cable ends.
3. Turn on the calculator, and start the DATAMATE program. Press **CLEAR** to reset the program.
4. Set up the calculator and interface for the temperature probe.
 - a. Select SETUP from the main screen.
 - b. If the calculator displays a temperature probe in CH 1, proceed directly to Step 5. If it does not, continue with this step to set up your sensor manually.
 - c. Press **ENTER** to select CH 1.
 - d. Select TEMPERATURE from the SELECT SENSOR menu.
 - e. Select the temperature probe you are using (in °C) from the TEMPERATURE menu.

Energy Content of Foods *continued*

5. Set up the data-collection mode.
 - a. To select MODE, press \blacktriangle once and press ENTER .
 - b. Select TIME GRAPH from the SELECT MODE menu.
 - c. Select CHANGE TIME SETTINGS from the TIME GRAPH SETTINGS menu.
 - d. Enter “6” as the time between samples in seconds.
 - e. Enter “100” as the number of samples. The length of the data collection will be 10 minutes.
 - f. Select OK to return to the setup screen.
 - g. Select OK again to return to the main screen.
6. Obtain a piece of one of the two foods assigned to you and a food holder like the one shown in **Figure 1**. Find and record the initial mass of the food sample and food holder. **CAUTION:** *Do not eat or drink in the laboratory.*
7. Determine and record the mass of an empty can. Obtain cold water from your teacher, and add 50 mL of it to the can. Determine and record the mass of the can and water.
8. Set up the apparatus as shown in **Figure 1**. Use a ring and stirring rod to suspend the can about 2.5 cm (1 in.) above the food sample. Use a utility clamp to suspend the temperature probe in the water. The probe should not touch the bottom of the can. Remember that the temperature probe must be in the water for at least 30 seconds before you complete Step 9.

DATA TABLE 1

Food sample 1:			
Initial mass of food sample and holder:			
Mass of empty can:		Mass of can and water:	
Food sample 2:			
Initial mass of food sample and holder:			
Mass of empty can:		Mass of can and water:	

DATA COLLECTION

9. Select START to begin collecting data. Record the initial temperature of the water, T_1 , in Data Table 2 (round to the nearest 0.1°C). **Note:** You can monitor temperature in the upper-right corner of the real-time graph displayed on the calculator screen. Remove the food sample from under the can, and use a wooden splint to light it. Quickly place the burning food sample directly under the center of the can. Allow the water to be heated until the food sample stops burning.
10. Continue stirring the water until the temperature stops rising. Record this maximum temperature, T_2 . Data collection will stop after 10 minutes (or press the $\text{STO} \blacktriangleright$ key to stop *before* 10 minutes have elapsed).
11. Determine and record the final mass of the food sample and food holder.

Energy Content of Foods *continued*

12. To confirm the initial (T_1) and final (T_2) values you recorded earlier, examine the data points along the curve on the displayed graph. As you move the cursor right or left, the time (X) and temperature (Y) values of each data point are displayed below the graph.
13. Press **ENTER** to return to the main screen. Select START to repeat the data collection for the second food sample. Use a new 50 mL portion of cold water. Repeat Steps 6–12.
14. When you are done, place burned food, used matches, and partially burned wooden splints in the container provided by the teacher.

DATA TABLE 2

Food sample 1:			
T_1 :		T_2 :	Final mass of sample and holder:
Food sample 2:			
T_1 :		T_2 :	Final mass of sample and holder:

Analysis

1. **Organizing data** Find the mass of water heated for each sample. _____

2. **Organizing data** Find the change in temperature of the water, ΔT , for each sample. _____

3. **Organizing data** Find the mass (in g) of each food sample burned. _____

4. **Analyzing Results** Calculate the heat absorbed by the water, q , using the equation

$$q = C_p m \Delta T$$

where q is heat, C_p is the specific heat, m is the mass of water, and ΔT is the change in temperature. For water, C_p is 4.18 J/g°C. Convert your final answer to units of kJ. _____

Energy Content of Foods *continued*

5. Analyzing Results Use the results of the previous two steps to calculate the energy content (in kJ/g) of each food sample. _____

DATA TABLE 3

Food sample 1:			
Mass of water heated:	g	Temperature change, ΔT :	$^{\circ}\text{C}$
Mass of food burned:	g	Heat, q :	kJ
<i>Energy content of food sample:</i>			kJ/g
Food sample 2:			
Mass of water heated:	g	Temperature change, ΔT :	$^{\circ}\text{C}$
Mass of food burned:	g	Heat, q :	kJ
<i>Energy content of food sample:</i>			kJ/g

Conclusions

1. Evaluating results Record your results and the results of other groups in the Class Results Table below. Which food had the highest energy content? Which had the lowest energy content? _____

CLASS RESULTS TABLE

Marshmallows	Peanuts	Cashews	Popcorn
kJ/g	kJ/g	kJ/g	kJ/g
kJ/g	kJ/g	kJ/g	kJ/g
kJ/g	kJ/g	kJ/g	kJ/g
kJ/g	kJ/g	kJ/g	kJ/g
kJ/g	kJ/g	kJ/g	kJ/g

Average for each food type:

kJ/g	kJ/g	kJ/g	kJ/g
------	------	------	------

2. Evaluating results Food energy is often expressed in a unit called a Calorie (or dietary calorie). There are 4.18 kJ in one Calorie. Based on the class average for popcorn, calculate the number of Calories in a 50.0 g package of popcorn. _____

Energy Content of Foods *continued*

3. Evaluating results Two of the foods in the experiment have a high fat content (peanuts and cashews), and two have a high carbohydrate content (marshmallows and popcorn). From your results, what generalization can you make about the relative energy content of fats and carbohydrates? _____

4. Evaluating results Based on the data you and your classmates collected, which of the four foods tested would you suggest to send on the NASA space mission?

Extensions

1. Applying results If you were packing for a mountain hike, what kind of snacks would you bring along? Why? _____

2. Critiquing methods Was all of the heat given off by the burning food sample transferred to the water in the can? How could this experiment be improved to account for all the heat given off when the food sample was burned?

3. Applying results Listed on the following page are possible nutrition labels for each of the food samples that you tested. Based on the data you and your classmates obtained in this lab, determine which of these labels is accurate and which is not. If you find a label to be incorrect, explain your reasoning.

Energy Content of Foods *continued*

MARSHMALLOWS

Nutrition Facts	
Serving Size 1 ounce	
Servings Per Container 6	
Amount per serving	
Calories 260	Calories from Fat 160
% Daily Value	
Total Fat 18g	13%
Saturated Fat 5g	27%
Cholesterol 0mg	0%
Sodium 260mg	11%
Total Carbohydrate 23g	8%
Dietary Fiber 1g	11%
Sugars 18g	
Protein 1g	

PEANUTS

Nutrition Facts	
Serving Size 1 ounce	
Servings Per Container 16	
Amount per serving	
Calories 165	Calories from Fat 125
% Daily Value	
Total Fat 14g	20%
Saturated Fat 1.9g	10%
Cholesterol 0mg	0%
Sodium 122mg	5%
Total Carbohydrate 5g	2%
Dietary Fiber 1g	4%
Sugars 2g	
Protein 8g	

CASHEWS

Nutrition Facts	
Serving Size 1 ounce	
Servings Per Container 16	
Amount per serving	
Calories 80	Calories from Fat 26
% Daily Value	
Total Fat 3g	4%
Saturated Fat 0.5g	3%
Cholesterol 0mg	0%
Sodium 177mg	7%
Total Carbohydrate 8g	3%
Dietary Fiber 2g	8%
Sugars 2g	
Protein 5g	

POPCORN

Nutrition Facts	
Serving Size 1 cup	
Servings Per Container 8	
Amount per serving	
Calories 30	Calories from Fat 0
% Daily Value	
Total Fat 0.3g	*0%
Saturated Fat 0g	*0%
Cholesterol 0mg	0%
Sodium 0mg	*0%
Total Carbohydrate 6g	2%
Dietary Fiber 1g	4%
Sugars 2g	
Protein 0g	
*Less than 1% of US RDA	