

TEACHER INFORMATION

Mapping the Ocean Floor

1. The student pages with complete instructions for data collection using LabQuest App, Logger *Pro* (computers), EasyData or DataMate (calculators), and DataPro (Palm handhelds) can be found on the CD that accompanies this book. See *Appendix A* for more information.
2. The boxes should be at least 50 cm wide.
3. Best results are obtained when the boxes are located touching each other or at least 50 cm apart.
4. The tabletops should be smooth and flat.
5. Boards longer than one meter can be used. They should provide rigid support for the Motion Detector.
6. The hidden arrangement can be located behind a counter or demonstration table.

ANSWERS TO QUESTIONS

For Sample Answers to the questions in this lab, please contact Vernier Software and Technology at swnanswers@vernier.com

ACKNOWLEDGEMENT

We wish to thank Don Volz and Sandy Sapatka for their help in developing and testing this experiment.